

ProRen Bladder Management System

System Collection Unit

Only for use with the ProRen Bladder Management System

ProRen FLO™

Developed for use in **Zero Gravity** for **Astronauts**.

Omni's patented, non-invasive, sensor-driven technology actively pulls urine away from the body, keeping the user dry, providing advanced incontinence and wound care management



ProRen FLO™ stays comfortably in place for ambulatory or non-ambulatory users day and night.

Quick Disconnect

Shown with System Control Unit, Waist Collection Chamber and System Chamber Unit.

Multiple collection chamber models are available:

- Waist Chamber
- Thigh Chamber
- Leg Chamber
- Bedside Chamber



Omni+

Designed for males with urinary incontinence



ProRen Bladder Management System

System Collection Unit

Only for use with the ProRen Bladder Management System

ProRen FLO™

The ProRen FLO™ solves the problems associated with traditional external catheters

The ProRen FLO™ uses an odor proof, hydro-block filter that enables air to freely flow into the collection unit. Urine sensors **automatically** activate the control unit, pulling urine completely away from the skin and **keeping the user dry**.

By allowing the collection unit to stay **fully drained**, the risk of developing urinary tract infections, skin macerations, penile erosion, and other possibly life threatening conditions is **drastically reduced**.

By keeping the skin dry, the collection unit remains **secure**.

The System Chamber Unit (hose/sensor) **easily detaches**, allowing the use of intermittent catheters.

With **Standard External Catheters**, there is no way for air to flow into the catheter, creating a vacuum which *prevents the urine from draining*, **trapping** the urine in the external catheter **against the skin**.

This can result in developing life threatening conditions such as **skin maceration**, **urinary tract infections**, and **penile erosion**. It also causes **accidental detachment**.

The ProRen FLO™ comes in a variety of configurations and sizes, ensuring men, including those who were previously hard to fit, will find an option that fits well and stays securely in place.

The ProRen FLO™ is only for use with the ProRen Bladder Management System.

